
Product Specification

Getinge Instrument Tray SPRI SPRI I, SPRI II, SPRI III Accessory

The Instrument Tray SPRI is designed for instrument circulation within a Central Sterile Services Department (CSSD) and Theatre Sterile Supplies Unit (TSSU), as well as at sub-sterilization units.

The instrument trays are intended for use through all phases of the instrument processing. The trays are are stackable, easy to clean and disinfect in a washer-disinfector. They are designed for fast self-drying. The trays are made of stainless steel electro-polished wire mesh with gaps of 6.5 mm (0.26"); large enough to allow the flow of air and water, yet small enough to prevent instrument damage.



The product range consists of three standard sizes, all adapted to fit inside modular wire baskets of the SPRI and other ISO standard baskets. Additionally, specialized instrument trays can be combined with a basic instrument tray set-up, covering all instruments used during any surgery. The trays are jig-welded to maintain size and shape even during heavyduty work.

Key features

- Stackable
- · Electrolytic polish for a smooth, clean surface
- Fulfills the hygiene demands of the health care sector and is autoclavable

Specifications

| | Instrument Tray SPRI I | Instrument Tray SPRI II | Instrument Tray SPRI III |
|--------------------------|---|---|---|
| Article No. | 6003000701 | 6003000702 | 6003000703 |
| Size L x W x H mm (") | 450 x 340 x 70 (57.09 x 13.39 x 2.76) | 340 x 250 x 70 (13.39 x 9.84 x 2.76) | 250 x 170 x 70 (9.84 x 6.69 x 2.76) |
| Weight kg (lb) | 1.14 (2.51) | 0.74 (1.63) | 0.4 (0.88) |
| Frame and mesh | Stainless steel AISI.304 Electro-polished wire | Stainless steel AISI.304 Electro-polished wire | Stainless steel AISI.304 Electro-polished wire |
| Mesh size mm (") | 6.5 (0.26) | 6.5 (0.26) | 6.5 (0.26) |

Standards and codes


This product is classified as a Class I medical device and complies with the requirements in the applicable standards and directives.



Accessories for trays

Lids

The lids are designed as a wire net surrounded by a rigid wire frame. The lid is attached with hooks at the top of the tray, and locked down with a spring wire. There are three lid sizes, one fitting each tray.



Specifications

| | Lid for Instrument Tray SPRI I | Lid for Instrument Tray SPRI II | Lid for Instrument Tray SPRI III |
|----------------------|--|--|--|
| Article No. | 477345200 | 477345300 | 498637400 |
| Size L x W mm (") | 450 x 340 (17.72 x 13.39) | 340 x 250 (13.39 x 9.84) | 250 x 170 (9.84 x 6.69) |
| Weight kg (lb) | 0.64 (1.41) | 0.38 (0.84) | 0.22 (0.49) |
| Frame and mesh | Stainless steel Electro-polished wire | Stainless steel Electro-polished wire | Stainless steel Electro-polished wire |

Distributed by: Getinge Disinfection AB · Ljungadalsgatan 11 · 352 46 Växjö · SWEDEN · +46 10 335 98 00